

General Purpose Water Bath, 10 Liter

WBE10A1

Key Specifications

Working Temperature:	Ambient +5° to 99°C
Temperature Stability:	±0.1°C
Power Requirements:	120V, 60 Hz or 220-240V, 50/60 Hz
Reservoir Capacity (gallons/liters):	2.64(10)
Working Access (L x W x H):	10.6 x 11.6 x 6" (26.9 x 29.5 x 15.2 cm)
Overall dimensions (L x W x H):	15.5 x 17 x 12" (39.3 x 43.2 x 30.5 cm)
Heater Wattage:	1000 W



Large, easy to read display provides an enhanced user experience and the ability to monitor temperature at a distance

Key Features

- Calibration offset feature
- Color TFT (thin-film transistor LCD) Display with adjustable brightness
- Five (5) Preset Programs allow the user to set a temperature and duration
- Integrated independent timer
- Clear, hinged high-gable cover
- Redundant safety thermostat
- Recessed handles

 **PolyScience®**
Temperature Control Solutions

PolyScience • 6600 W. Touhy Ave. • Niles, IL 60714 • www.PolyScience.com
Toll-Free +1.800.229.7569 • Phone +1.847.647.0611 • sales@PolyScience.com